

Product datasheet for **GA111127**

Human PCDHAC1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDHAC1 CRISPRa kit - CRISPR gene activation of human protocadherin alpha subfamily C, 1
Symbol:	PCDHAC1
Components:	GA111127G1, PCDHAC1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111127G2, PCDHAC1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111127G3, PCDHAC1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H158
Synonyms:	PCDH-ALPHA-C1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)