

Product datasheet for **GA111133**

Human PCDHA7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDHA7 CRISPRa kit - CRISPR gene activation of human protocadherin alpha 7
Symbol:	PCDHA7
Components:	GA111133G1, PCDHA7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111133G2, PCDHA7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111133G3, PCDHA7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UN72</u>
Synonyms:	CNR4; CNRN4; CNRS4; CRNR4; PCDH-ALPHA7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)