

Product datasheet for **GA111137**

Human PCDHA3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDHA3 CRISPRa kit - CRISPR gene activation of human protocadherin alpha 3
Symbol:	PCDHA3
UniProt ID:	Q9Y5H8
Synonyms:	PCDH-ALPHA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)