

Product datasheet for **GA111155**

Human Leiomodlin 3 (LMOD3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LMOD3 CRISPRa kit - CRISPR gene activation of human leiomodlin 3
Symbol:	LMOD3
Components:	GA111155G1, Leiomodlin 3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111155G2, Leiomodlin 3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111155G3, Leiomodlin 3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q0VAK6</u>
Synonyms:	NEM10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)