

Product datasheet for **GA111167**

Human SERTAD4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SERTAD4 CRISPRa kit - CRISPR gene activation of human SERTA domain containing 4
Symbol:	SERTAD4
Components:	GA111167G1, SERTAD4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111167G2, SERTAD4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111167G3, SERTAD4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NUC0</u>
Synonyms:	DJ667H12.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)