


Product datasheet for **GA111175**

Human WDR45L (WDR45B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR45B CRISPRa kit - CRISPR gene activation of human WD repeat domain 45B
Symbol:	WDR45B
Components:	GA111175G1, WDR45L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111175G2, WDR45L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111175G3, WDR45L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5MNZ6
Synonyms:	WDR45L; WIPI-3; WIPI3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)