

Product datasheet for **GA111178**

Human PARD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PARD3 CRISPRa kit - CRISPR gene activation of human par-3 family cell polarity regulator
Symbol:	PARD3
UniProt ID:	<u>Q8TEW0</u>
Synonyms:	ASIP; Baz; PAR3; PAR3alpha; PARD-3; PARD3A; PPP1R118; SE2-5L16; SE2-5LT1; SE2-5T2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)