

Product datasheet for **GA111187**

Human Leukotriene B4 Receptor 2 (LTB4R2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LTB4R2 CRISPRa kit - CRISPR gene activation of human leukotriene B4 receptor 2
Symbol:	LTB4R2
UniProt ID:	<u>Q9NPC1</u>
Synonyms:	BLT2; BLTR2; JULF2; KPG_004; LTB4-R 2; LTB4-R2; NOP9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)