

Product datasheet for **GA111189**

Human RPRM activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RPRM CRISPRa kit - CRISPR gene activation of human reprimo, TP53 dependent G2 arrest mediator homolog
Symbol:	RPRM
UniProt ID:	Q9NS64
Synonyms:	REPRIMO
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)