

Product datasheet for **GA111192**

Human **KCNQ5** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNQ5 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily Q member 5
Symbol:	KCNQ5
UniProt ID:	Q9NR82
Synonyms:	Kv7.5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)