

Product datasheet for **GA111204**

Human TMPRSS4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMPRSS4 CRISPRa kit - CRISPR gene activation of human transmembrane serine protease 4
Symbol:	TMPRSS4
UniProt ID:	Q9NRS4
Synonyms:	CAPH2; MT-SP2; TMPRSS3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)