

Product datasheet for **GA111206**

Human Twinkle (TWNK) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TWNK CRISPRa kit - CRISPR gene activation of human twinkle mtDNA helicase
Symbol:	TWNK
Components:	GA111206G1, Twinkle gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111206G2, Twinkle gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111206G3, Twinkle gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96RR1</u>
Synonyms:	ATXN8; IOSCA; MTDP57; PEO; PEO1; PEOA3; PRLTS5; SANDO; SCA8; TWINL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)