

Product datasheet for **GA111212**

Human **KCNK12** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNK12 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 12
Symbol:	KCNK12
UniProt ID:	Q9HB15
Synonyms:	K2p12.1; THIK-2; THIK2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)