

Product datasheet for **GA111219**

Human NRIP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NRIP3 CRISPRa kit - CRISPR gene activation of human nuclear receptor interacting protein 3
Symbol:	NRIP3
UniProt ID:	Q9NQ35
Synonyms:	C11orf14; NY-SAR-105
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)