

Product datasheet for **GA111225**

Human **SLC2A4RG** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC2A4RG CRISPRa kit - CRISPR gene activation of human SLC2A4 regulator
Symbol:	SLC2A4RG
UniProt ID:	Q9NR83
Synonyms:	GEF; HDBP-1; HDBP1; Si-1-2; Si-1-2-19
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)