

## Product datasheet for **GA111242**

### Human TM9SF3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TM9SF3 CRISPRa kit - CRISPR gene activation of human transmembrane 9 superfamily member 3
Symbol:	TM9SF3
UniProt ID:	<a href="#">Q9HD45</a>
Synonyms:	EP70-P-iso; SMBP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)