

Product datasheet for GA111250

Human WHIP (WRNIP1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WRNIP1 CRISPRa kit - CRISPR gene activation of human WRN helicase interacting protein 1
Symbol:	WRNIP1
UniProt ID:	Q96S55
Synonyms:	bA420G6.2; WHIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)