

Product datasheet for **GA111267**

Human SMARCAD1 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | SMARCAD1 CRISPRa kit - CRISPR gene activation of human SWI/SNF-related, matrix-associated actin-dependent regulator of chromatin, subfamily a, containing DEAD/H box 1 |
| Symbol: | SMARCAD1 |
| UniProt ID: | Q9H4L7 |
| Synonyms: | ADERM; ETL1; HEL1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)