


Product datasheet for **GA111270**

Human DHX33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DHX33 CRISPRa kit - CRISPR gene activation of human DEAH-box helicase 33
Symbol:	DHX33
Components:	GA111270G1, DHX33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111270G2, DHX33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111270G3, DHX33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H6R0</u>
Synonyms:	DDX33
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)