

Product datasheet for **GA111277**

Human GPR108 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GPR108 CRISPRa kit - CRISPR gene activation of human G protein-coupled receptor 108
Symbol:	GPR108
Components:	GA111277G1, GPR108 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111277G2, GPR108 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111277G3, GPR108 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NPR9</u>
Synonyms:	LISTR2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)