

Product datasheet for **GA111280**

Human **DUS3L** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	DUS3L CRISPRa kit - CRISPR gene activation of human dihydrouridine synthase 3 like
Symbol:	DUS3L
UniProt ID:	Q96G46
Synonyms:	DUS3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)