

## Product datasheet for **GA111293**

### Human MFF activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                      |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                         |
| Description:  | MFF CRISPRa kit - CRISPR gene activation of human mitochondrial fission factor                                                                                                                                                                                                                       |
| Symbol:       | MFF                                                                                                                                                                                                                                                                                                  |
| Components:   | GA111293G1, MFF gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111293G2, MFF gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111293G3, MFF gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9GZY8</a></u>                                                                                                                                                                                                                                                                        |
| Synonyms:     | C2orf33; EMPF2; GL004                                                                                                                                                                                                                                                                                |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                              |

#### Product images:



[View online »](#)