


Product datasheet for **GA111304**

Human SHD activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SHD CRISPRa kit - CRISPR gene activation of human Src homology 2 domain containing transforming protein D
Symbol:	SHD
Components:	GA111304G1, SHD gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111304G2, SHD gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111304G3, SHD gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96IW2</u>
Synonyms:	src homology 2 domain-containing transforming protein D; Src homology 2 domain containing transforming protein D
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)