

Product datasheet for **GA111318**

Human PSMG2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMG2 CRISPRa kit - CRISPR gene activation of human proteasome assembly chaperone 2
Symbol:	PSMG2
Components:	GA111318G1, PSMG2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111318G2, PSMG2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111318G3, PSMG2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q969U7</u>
Synonyms:	CLAST3; HCCA3; HsT1707; MDS003; PAC2; TNFSF5IP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)