

Product datasheet for **GA111334**

Human ACKR3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ACKR3 CRISPRa kit - CRISPR gene activation of human atypical chemokine receptor 3
Symbol:	ACKR3
UniProt ID:	P25106
Synonyms:	CMKOR1; CXC-R7; CXCR-7; CXCR7; GPR159; RDC-1; RDC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)