


Product datasheet for **GA111340**

Human VPS35L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	VPS35L CRISPRa kit - CRISPR gene activation of human VPS35 endosomal protein sorting factor like
Symbol:	VPS35L
Components:	GA111340G1, VPS35L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111340G2, VPS35L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111340G3, VPS35L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z3J2</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)