

Product datasheet for **GA111346**

Human TWSG1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TWSG1 CRISPRa kit - CRISPR gene activation of human twisted gastrulation BMP signaling modulator 1
Symbol:	TWSG1
UniProt ID:	Q9GZX9
Synonyms:	TSG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)