

Product datasheet for **GA111355**

Human MCG10 (PCBP4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCBP4 CRISPRa kit - CRISPR gene activation of human poly(rC) binding protein 4
Symbol:	PCBP4
UniProt ID:	P57723
Synonyms:	CBP; LIP4; MCG10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)