

Product datasheet for **GA111393**

Human LYRM4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LYRM4 CRISPRa kit - CRISPR gene activation of human LYR motif containing 4
Symbol:	LYRM4
Components:	GA111393G1, LYRM4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111393G2, LYRM4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111393G3, LYRM4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9HD34
Synonyms:	C6orf149; CGI-203; COXPD19; ISD11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)