


Product datasheet for **GA111394**

Human MRPL47 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MRPL47 CRISPRa kit - CRISPR gene activation of human mitochondrial ribosomal protein L47
Symbol:	MRPL47
Components:	GA111394G1, MRPL47 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111394G2, MRPL47 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111394G3, MRPL47 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9HD33
Synonyms:	CGI-204; L47mt; MRP-L47; NCM1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)