

Product datasheet for **GA111418**

Human Pellino 2 (PELI2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PELI2 CRISPRa kit - CRISPR gene activation of human pellino E3 ubiquitin protein ligase family member 2
Symbol:	PELI2
UniProt ID:	Q9HAT8
Synonyms:	pellino 2; pellino homolog 2 (Drosophila)
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)