

Product datasheet for **GA111419**

Human Pellino 1 (PELI1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PELI1 CRISPRa kit - CRISPR gene activation of human pellino E3 ubiquitin protein ligase 1
Symbol:	PELI1
UniProt ID:	Q96FA3
Synonyms:	DKFZp686C18116; MGC50990
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)