


Product datasheet for **GA111434**

Human NIPAL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NIPAL3 CRISPRa kit - CRISPR gene activation of human NIPA like domain containing 3
Symbol:	NIPAL3
Components:	GA111434G1, NIPAL3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111434G2, NIPAL3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111434G3, NIPAL3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P499</u>
Synonyms:	DJ462O23.2; NPAL3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)