

## Product datasheet for **GA111446**

### Human SLC45A4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SLC45A4 CRISPRa kit - CRISPR gene activation of human solute carrier family 45 member 4
Symbol:	SLC45A4
Components:	GA111446G1, SLC45A4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111446G2, SLC45A4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111446G3, SLC45A4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q5BKX6</a></u>
Synonyms:	KIAA1126
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)