

Product datasheet for **GA111453**

Human ARFGEF3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARFGEF3 CRISPRa kit – CRISPR gene activation of human ARFGEF family member 3
Symbol:	ARFGEF3
Components:	GA111453G1, ARFGEF3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111453G2, ARFGEF3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111453G3, ARFGEF3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5TH69
Synonyms:	A7322; BIG3; C6orf92; dJ171N11.1; PPP1R33
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

