

Product datasheet for GA111462

Human KIR2DL5A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIR2DL5A CRISPRa kit - CRISPR gene activation of human killer cell immunoglobulin like receptor, two Ig domains and long cytoplasmic tail 5A
Symbol:	KIR2DL5A
UniProt ID:	Q8N109
Synonyms:	CD158F; KIR2DL5; KIR2DL5.1; KIR2DL5.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)