

Product datasheet for **GA111464**

Human PBXIP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PBXIP1 CRISPRa kit - CRISPR gene activation of human PBX homeobox interacting protein 1
Symbol:	PBXIP1
UniProt ID:	<u>Q96AQ6</u>
Synonyms:	HPIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)