


Product datasheet for **GA111477**

Human CLTRN activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CLTRN CRISPRa kit - CRISPR gene activation of human collectrin, amino acid transport regulator
Symbol:	CLTRN
Components:	GA111477G1, CLTRN gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111477G2, CLTRN gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111477G3, CLTRN gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9HBJ8
Synonyms:	NX-17; NX17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)