

Product datasheet for **GA111489**

Human TMIGD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMIGD3 CRISPRa kit - CRISPR gene activation of human transmembrane and immunoglobulin domain containing 3
Symbol:	TMIGD3
Components:	GA111489G1, TMIGD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111489G2, TMIGD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111489G3, TMIGD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P33765</u>
Synonyms:	AD026
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

