

Product datasheet for **GA111498**

Human TENM2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TENM2 CRISPRa kit - CRISPR gene activation of human teneurin transmembrane protein 2
Symbol:	TENM2
Components:	GA111498G1, TENM2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111498G2, TENM2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111498G3, TENM2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	ODZ2; ten-2; TEN-M2; TNM2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)