


Product datasheet for **GA111503**

Human TMCC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMCC3 CRISPRa kit - CRISPR gene activation of human transmembrane and coiled-coil domain family 3
Symbol:	TMCC3
Components:	GA111503G1, TMCC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111503G2, TMCC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111503G3, TMCC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9ULS5</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)