

Product datasheet for **GA111515**

Human **LRRC47** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC47 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 47
Symbol:	LRRC47
Components:	GA111515G1, LRRC47 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111515G2, LRRC47 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111515G3, LRRC47 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N1G4</u>
Synonyms:	RP1-286D6.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)