

Product datasheet for **GA111523**

Human USP31 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	USP31 CRISPRa kit - CRISPR gene activation of human ubiquitin specific peptidase 31
Symbol:	USP31
Components:	GA111523G1, USP31 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111523G2, USP31 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111523G3, USP31 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q70CQ4</u>
Synonyms:	RP11-20G6.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)