

Product datasheet for **GA111525**

Human PLEKHG1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHG1 CRISPRa kit - CRISPR gene activation of human pleckstrin homology and RhoGEF domain containing G1
Symbol:	PLEKHG1
Components:	GA111525G1, PLEKHG1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111525G2, PLEKHG1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111525G3, PLEKHG1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9ULL1</u>
Synonyms:	ARHGEF41
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)