

## Product datasheet for **GA111536**

### Human KIAA1239 (NWD2) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | NWD2 CRISPRa kit - CRISPR gene activation of human NACHT and WD repeat domain containing 2                                                                                                                                                                                                                          |
| Symbol:       | NWD2                                                                                                                                                                                                                                                                                                                |
| Components:   | GA111536G1, KIAA1239 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111536G2, KIAA1239 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111536G3, KIAA1239 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q9ULI1</a>                                                                                                                                                                                                                                                                                              |
| Synonyms:     | KIAA1239                                                                                                                                                                                                                                                                                                            |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)