

Product datasheet for **GA111538**

Human LRFN2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRFN2 CRISPRa kit - CRISPR gene activation of human leucine rich repeat and fibronectin type III domain containing 2
Symbol:	LRFN2
Components:	GA111538G1, LRFN2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111538G2, LRFN2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111538G3, LRFN2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9ULH4
Synonyms:	FIGLER2; KIAA1246; SALM1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)