


## Product datasheet for **GA111549**

### Human COG6 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | COG6 CRISPRa kit - CRISPR gene activation of human component of oligomeric golgi complex 6                                                                                                                                                                                                              |
| Symbol:       | COG6                                                                                                                                                                                                                                                                                                    |
| Components:   | GA111549G1, COG6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111549G2, COG6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111549G3, COG6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9Y2V7</a></u>                                                                                                                                                                                                                                                                           |
| Synonyms:     | CDG2L; COD2; SHNS                                                                                                                                                                                                                                                                                       |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)