

## Product datasheet for **GA111550**

### Human GPR158 activation kit by CRISPRa

#### Product data:

|               |                                                                                     |
|---------------|-------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                        |
| Description:  | GPR158 CRISPRa kit - CRISPR gene activation of human G protein-coupled receptor 158 |
| Symbol:       | GPR158                                                                              |
| UniProt ID:   | <u><a href="#">Q5T848</a></u>                                                       |
| Synonyms:     | RP11-59G22.1                                                                        |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)             |

#### Product images:



[View online >](#)