


Product datasheet for **GA111572**

Human WDR35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR35 CRISPRa kit - CRISPR gene activation of human WD repeat domain 35
Symbol:	WDR35
Components:	GA111572G1, WDR35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111572G2, WDR35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111572G3, WDR35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9P2L0</u>
Synonyms:	CED2; IFT121; IFTA1; SRTD7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)