


Product datasheet for **GA111573**

Human PTCHD2 (DISP3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DISP3 CRISPRa kit - CRISPR gene activation of human dispatched RND transporter family member 3
Symbol:	DISP3
Components:	GA111573G1, PTCHD2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111573G2, PTCHD2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111573G3, PTCHD2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9P2K9
Synonyms:	PTCHD2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)